

## OXYGEN FLASK CALIBRATION PROCEDURE

### **CLEANING**

Wash each flask and stopper thoroughly inside and outside with tap water x3 and then distilled water x3. Once clean, **DO NOT TOUCH THE FLASK WITH YOUR FINGERS – WEAR GLOVES.**

### **WEIGHING**

#### **A. Pre-weighing procedure**

- Thermometers need to be within a ½ degree of each other.
- Ensure balance is on and stable.
  - The balance will need to have a range up to 300 grams and weigh to the thousandths place (1 milligram).
  - Ideally, each flask can be known to  $\pm 0.003$  mL.
  - Use the METTLER TOLEDO CO2 lab's analytical balance.

#### **B. Dry weighing (g)**

- While wearing gloves, place the flask with its stopper on onto the weigh pan and weigh. Allow the scale to stabilize and record the weight.
  - Weigh at least twice to confirm weight.
  - Record dry weight

#### **C. Wet weighing (g)**

- Have distilled water that has been sitting out for at least two hours.
- Rinse the stopper and the neck of the flask.
- Pinch noodle, place at bottom of flask at an angle and slowly fill the flask with distilled water (overflow).
- Pinch tube to pull noodle out.
- Place the stopper onto the flask. Ensure there are no trapped bubbles in the flask.
- Pour off any excess water that was displaced by the stopper. Carefully, wipe excess water off the flask and around the stopper.
  - Use spatula and Kimwipes to get down around stopper in neck.
- While wearing gloves, place the flask and its stopper on the weigh pan and weigh. Allow the scale to stabilize and record the weight.
  - Weigh at least twice to confirm weight.
  - Record wet weight.

#### **D. Temperature (Celsius)**

- Remove the flask from the weigh pan. Insert a thermometer into the flask. Record the temperature to the tenths place (i.e. 22.7) after two minutes.
  - Can also have a thermometer in a beaker of the same DI water next to setup and use that temperature.
  - Record temperature.

## CALCULATIONS

**A.  $rw = 0.999842594 + 0.00006793952 * (\text{temperature}) - 0.00000909529 * (\text{temperature})^2 + 0.0000001001685 * (\text{temperature})^3 - 0.000000001120083 * (\text{temperature})^4 + 0.000000000006536332 * (\text{temperature})^5$**

**B.  $Vol (@RT) = (\text{buoyancy correction}) * ((\text{filled weight}) - (\text{dry weight})) / (rw)$**   
buoyancy correction = 1.00105

**C. Vol (20C) = (Vol (@RT))\*(1 + 0.00001\*(20 – (temperature)))**